



Polishing single-column cartridge, part no. 10038
Deionization single-column cartridge, part no. 10318

MAIN FEATURES

- Designed for Compact UP ultrapure water systems
- Tool-free installation
- Durable PP housing
- High-quality ion exchange resins
- Compact dimensions: W 8.5 × D 8.0 × H 37.0 cm

SPECIFICATION FOR SINGLE-COLUMN CARTRIDGES

Part No.	Resource*	Replacement interval	Feed water**	Corresponding systems
10318	1500l	1 year	≤800 µS/cm / RO TDS ≥98%	Compact Trace/ HPLC/ Bio (produced from July 2026)
10038	3000 l	1 year	≤1 µS/cm	Compact Trace/ HPLC/ Bio (produced from July 2026)

* Depends on the feed water ionic content and the concentration of dissolved CO₂.

** Feed water properties for optimal resource.

Storage period not longer than 2 years.

Store in sealed original packaging in a dry, ventilated area at 5–30 °C.